

Shoprider

PUMA

ASK FOR THE BRAND SHOPRIDER



PUMA 10



PUMA 14HD

SHOPRIDER[®]
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®

www.shoprider.com.au

mobility

PUMA

The SHOPRIDER™ PUMA has been ergonomically engineered to provide a robust combination of style and comfort. The PUMA is fitted with powerful smooth motors, mid wheel manoeuvrability, and comfortable seating. The smooth operation of the front and rear elastometers combined with mid wheel drive means easy operation on a variety of terrains and hard surfaces. It's ideal for both inside and outside use on even ground and can travel long distances between recharging.

SPECIFICATIONS

	PUMA10	PUMA14HD
Length	1060mm	1060mm
Width	620mm	640mm
Height	1160mm	1180mm
Battery Capacity	2 x 36Amp	2 x 75Amp
#Maximum Range	up to 20km	up to 25km
Maximum Load Single User*	136Kgs	205Kgs
Charger on board	4Amp	8 Amp
Motor Size	0.8 Hp	1.34 Hp
Maximum Speed	8Kph	6Kph
*Safe Climbing Angle	8 Degrees	8 Degrees
*Max Climbing Angle	12 Degrees	9 Degrees
Turning Radius	570mm	600mm
Wheel Size:		
- Front Caster	2 x 175mm	2 x 175mm
- Middle (drive)	2 x 260mm	2 x 350mm
- Rear caster	2 x 175mm	2 x 175mm

FEATURES

Large Captain's Seat
 20" Seat PUMA 14HD
 Sleek body design
 Smooth operation
 Variable Seat height
 Controller interchangeable both arms
 Free wheel device
 Adjustable arm rests
 Option elevating leg rests
 Option of swing away leg rest
 Stamp supports
Comprehensive Warranty


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10kph on public roads or paths

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPRIDER®

www.shoprider.com.au

mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
 Or visit our Web Site at www.shoprider.com.au